USER MANUAL

RAPID TEST FOR THE DETECTION OF
BETA-LACTAM, TETRACYCLINE, CHLORAMPHENICOL
AND STREPTOMYCIN ANTIBIOTICS IN MILK

PRODUCT REFERENCE
KIT060

ENGLISH

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Version 9 - October 30, 2020
Designed and produced in Belgium
© Unisensor • All rights reserved


TABLE OF CONTENTS**ENGLISH**

HOW TO USE A MICROPIPETTE	3
SENSITIVITY TABLE A	4
CHAPTER 1 - INTRODUCTION	6
CHAPTER 2 - COMPOSITION OF THE KIT	7
CHAPTER 3 - GENERAL REMARKS	8
LIMITATION OF USE	8
CHAPTER 4 - OPERATING PROCEDURE	9
CHAPTER 5 - RESULTS	11
VISUAL INTERPRETATION	11
INSTRUMENTAL INTERRETATION: USING READSENSOR 2 / ANKAR-100	12
CHAPTER 6 - NEGATIVE & POSITIVE CONTROLS RECONSTITUTION	13

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ENGLISH

CHAPTER 1 - INTRODUCTION

4Sensor (KIT060) is a rapid qualitative test that allows you to simultaneously detect and discriminate the presence or absence of **β -lactams, Tetracyclines, (Dihydro)-Streptomycin and Chloramphenicol molecules** in a milk sample.

SUMMARY OF THE PROTOCOL

1. **Add 200 μ l of milk into one reagent microwell and mix to homogeneity**
2. **Incubate 5 minutes at 40°C**
3. **Dip one dipstick into each microwell**
4. **Continue incubating for 5 minutes at 40°C**
5. **Remove the dipstick and stop the reaction by removing the sample pad**
6. **Read the results**

REACTION MECHANISM

4Sensor (KIT060) is a competitive test involving specific receptors and generic monoclonal antibodies in one single operation.

The test requires the use of two components.

- The first component is a microwell containing predetermined amounts of both receptors and antibodies linked to gold particles.
- The second is a dipstick made up of a set of membranes with specific capture lines.

For a valid test, the upper red control line has to be visible after the second incubation.

The other four are specific "test" lines placed below the control line (see Figure A).

When the reagent from the microwell is reconstituted with a milk sample, both receptors and monoclonal antibodies will bind the corresponding analytes if present during the first 5 minutes incubation at 40°C.

Afterwards, when the dipstick is dipped into the milk, the liquid starts running vertically on the dipstick and passes through capture zones.

When the sample is free of antibiotics, a color development occurs at the specific capture lines, indicating the absence of the targeted analytes in milk sample. On the contrary, the presence of antibiotics in the sample will not cause the colored signal to appear at the specific capture lines.

CHAPTER 2 - COMPOSITION OF THE KIT

4Sensor Milk Kits contain everything needed to perform 96 measurements.

TEST TUBES

12 tubes each with 1 strip of 8 reagent microwells and 8 dipsticks.

POSITIVE CONTROLS

1 glass vial of POSITIVE CONTROL (**red cap**)


NEGATIVE CONTROLS

1 glass vial of NEGATIVE CONTROL (**green cap**)

MICROPIPETTE & TIPS

1 micropipette of 200 µl and disposable tips.

THE CERTIFICATE OF COMPLIANCE

USER MANUAL

ADDITIONAL MATERIAL NEEDED

- HeatSensor (40°C incubation, refer to the HeatSensor User Manual)
- ReadSensor 2 / ANKAR-100 (optional, refer to the ReadSensor 2 / ANKAR-100 User Manual)
- Distilled water

CHAPTER 3 - GENERAL REMARKS

- Store the kit at a temperature between 2°C and 8°C in a controlled fridge.
- Avoid repeated temperature variation.
- Avoid exposure of the reagents to moisture and light.
- Do not mix components of kits from different batches.
- Do not use the kit after the expiration date.
- Avoid touching the membrane on the dipstick.
- Avoid leak of reagent during sample addition and mixing step.
- Empty one tube before opening another tube.
- Close accurately the tube after use.
- The milk sample must be liquid and homogeneous. There can be neither clots or sedimentation phases.
- The ideal temperature of the milk sample is between 4 and 20°C.
- The temperature to perform the test is 40°C ± 1°C. Use the HeatSensor. Any other type of incubator is not appropriate to perform the 4Sensor assay. The use of other incubators that indicated in this document is the responsibility of the client and does not bind Unisensor (see HeatSensor user manual).
- The timing to perform the test is 5 minutes + 5 minutes. This timing must be respected, no extra time (second) is allowed.
- When drying, the color intensities of the lines will become sharper.
- When a positive result is recorded, the test result should be confirmed by performing two more tests.

LIMITATION OF USE

- The room temperature should be between 18°C and 24°C.
- Avoid to expose dipsticks to direct air flow like air conditioner, opened windows, etc. during migration.
- Respect recommended temperatures for milk samples storage (4°C to 8°C).
- Do not use other milks than fresh raw cow milk.
- Milk should be tested within 48 hours after milking.
- No frozen, transformed or powder milk can be tested.
- High fat samples might have an impact on the results.
- Level of somatic cells or bacterial presence might have an impact on results.
- To make **4Sensor (KIT060)** measurements reliable, pH of the milk must be between 6.5 and 7.2.
- Do not use in environment where relative humidity is outside the normal range (20%-80% at 20°C).
- Do not leave dipsticks outside the tube for a long period before use.

CHAPTER 4 - OPERATING PROCEDURE

PREPARATION

This procedure is described to easily run one single sample or a set of many samples. In that case, try to perform the test in cascade and avoid any delays when mixing reagent and milk but also when adding and removing dipsticks.

Make sure you have the same incubation time and temperature for each sample. You shouldn't test more than 8 samples at one time and, if there are more than 2 samples, you should use a multichannel pipette. With more than 8 samples we recommend to run series of maximum 8 samples.

1. Site preparation

Choose a clean and dry place to perform the test and wash and dry your hands before starting.

2. Start the incubator

Start the incubator (refer to the HeatSensor User Manual) and wait until the temperature has stabilized at 40°C.

3. Take the kit out of the fridge

Before opening the reagents, take the kit out of the fridge and wait until the temperature of the reagents reaches the ambient temperature. Meanwhile, read carefully the operating procedure. Determine how many samples are to be tested and write an identification number on each milk sample tube.

4. Open one plastic tube

Take out as many microwells as there are milk samples to be tested (the positive and negative controls included if necessary).


- If you do not intend to use all the 8 microwells, leave the sealing plastic foil on the unused ones.
- Put back unused microwells immediately into the white tube without damaging the dipsticks. Close and make sure it is tightly sealed.
- The plastic tubes should always be well closed after reagents have been taken out.

CHAPTER 4 - OPERATING PROCEDURE

5. Microwell(s) on incubator

Place the microwell(s) in the heating block which shows 40°C.

6. Transfer 200 µl of milk

Immediately transfer with a new tip 200 µl of milk into each of the microwells. Then mix by pipetting up and down 5 to 10 times.
Warning: when reagents and milk are in contact, the reaction begins.

Then, immediately push the start (run) button on the Heatsensor or use a timer. The 5-minutes countdown starts.

7. Place dipsticks into the microwells

During the 5 minutes incubation, open the same tube as before, take out as many dipsticks as there are analyses in progress and close the tube. Lay the dipsticks on a clean surface and identify them according to the milk sample IDs. When the 5 minutes are over, place manually the dipsticks into the corresponding microwells. The migration step of 5 minutes begins.

8. Sample pad removal

When the 5 minutes migration are over, take the dipsticks out of the microwells. Immediately remove the sample pad and proceed to the dipstick result interpretation.

NOTE

IF YOU ARE NOT PLANNING TO PERFORM ANY OTHER TEST WITHIN THE DAY, PUT EVERYTHING BACK INTO THE KIT AND STORE IT IN A FRIDGE AT A TEMPERATURE RANGING FROM 2°C TO 8°C.

CHAPTER 5 - RESULTS

VISUAL INTERPRETATION

1. Check whether the top control line is present. If it is not, consider the analysis as invalid and do not start (or continue) any interpretation.
2. When the top control line can be seen, interpret the four test lines as follows: examine one test line at a time and compare the intensity of the line color of the test line with the intensity of the line color of the control line (see Figure A). Start with the bottom line of β -lactam antibiotics for example.
 - If the test line is darker than the control line, the result is **NEGATIVE**, which means that, given the sensitivity of the test, the milk sample contains no antibiotics or antibiotics at a lower level than the value stated in the enclosed **TABLE A**.
 - If the test line is lighter or similar to the control line, the result is **POSITIVE**, which means that, given the sensitivity of the test, the milk sample contains antibiotics at or above the detection values stated in the enclosed **TABLE A**.
3. When you have interpreted one test line, do the same for the other lines.
4. If you hesitate, consider the sample as **POSITIVE**.
5. Write down your assessment on each of the dipsticks.
6. Dipsticks can be archived as a permanent record if required, by allowing dipstick to dry before storage.
7. Note: line color intensity will darken on drying.

CHAPTER 5 - RESULTS

INSTRUMENTAL INTERPRETATION: USING READSSENSOR 2 / ANKAR-100

1. Reading Method

Make sure you have entered the latest version available of the reading method into your reader (Check Unisensor Website for the latest appropriate version).

2. Reading

Refer to **ReadSensor 2 / ANKAR-100 - APP088** User Manual for detailed reading procedure.

- If the control is not visible (with or without the presence of test lines), do consider the test as not valid, don't interpret the result and directly run another test. The reader will confirm it by giving an INVALID result.
- You should read the dipstick within 5 minutes of performing the test.

On the basis of a ratio obtained between the test and control line color's intensities, the reader will interpret the result as follows:

RESULTS INTERPRETATION		
TEST RESULT	NEGATIVE (NEG)	POSITIVE (POS)
INTERPRETATION	Contamination absent or below the limit of detection	Contamination at or above the limit of detection

CHAPTER 6 - NEGATIVE & POSITIVE CONTROLS RECONSTITUTION


1. Take a glass vial of negative and positive controls

Take a glass vial of a negative (green cap) and a positive (red cap) control out of the bag.

Remove the cap and rubber stopper from the glass vials.

Avoid mixing or changing the respective stoppers of both positive and negative control vials.

2. Reconstitute with water

Add 1mL of distilled or deionized water at ambient temperature.

Replace the corresponding stopper on respective control vial and leave to stand for 1 minute and a half at ambient temperature.

3. Mix

Vigorously shake the vial during 10 seconds.

Check that the sample is properly dissolved. The material is now ready to use

Use the 200 μ L of control solution like any other milk samples for testing (see CHAPTER 4: Operating Procedure).

CHAPTER 7 - MILK POWDER DILUTION

In an appropriate bottle, mix 10 g of milk powder with 90 ml of warm (40°C) and distilled water.

For an optimal dilution, shake vigorously.

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DISCLAIMER

The information contained in this document is provided in good faith, and every reasonable effort is made to ensure that it is correct and up to date.

However, given the dynamic nature of our business, and changing rules and regulations, there may be delays, omissions or inaccuracies in information contained in this document. Accordingly, Unisensor does not warrant the accuracy and completeness of the information.

The information contained in this document may be changed or updated at any time without notice. The information provided was last updated as of October 30, 2020.

Any person relying on any of the information contained in this document or making any use of the information contained herein, shall do so at its own risk.

To the best of our knowledge and belief, the products mentioned herein do not infringe the intellectual property rights of any third party. The products may be covered by pending or issued patents, registered or unregistered trademarks, or similar intellectual property rights.

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